


Product Datasheet

Rat GLP-1 ELISA Kit (orb567419)

Description	Rat GLP-1 ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	GLP-1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P06883
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat GLP-1 (Glucagon Like Peptide 1) ELISA...