


Product Datasheet

Rat GLbeta ELISA Kit (orb567422)

Description	Rat GLbeta ELISA Kit
Reactivity	Rat
Range	0.156-10mIU/ml
Target	GL β
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	D3ZUM4
Sensitivity	0.094mIU/ml
Expiration Date	Please enquire.



Rat
GL β (Galactosidase,
Beta) ELISA Kit