


## Product Datasheet

### Rat HbA1C ELISA Kit (orb567443)

<b>Description</b>	Rat HbA1C ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	1.563-100ng/ml
<b>Target</b>	HbA1C
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P01946</b>
<b>Sensitivity</b>	0.938ng/ml
<b>Expiration Date</b>	Please enquire.



Rat  
HbA1C(Glycosylated  
Hemoglobin/Hemogl...