


Product Datasheet

Rat LH ELISA Kit (orb567514)

Description	Rat LH ELISA Kit
Reactivity	Rat
Range	0.313-20mIU/ml
Target	LH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01230
Sensitivity	0.188mIU/ml
Expiration Date	Please enquire.



Rat
LH(Luteinizing
Hormone)
ELISA Kit