


Product Datasheet

Rat CCL22 ELISA Kit (orb567526)

Description	Rat CCL22 ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	MDC/CCL22
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9QZU1
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
MDC/CCL22(Macrophage-Derived Chemoki...