


Product Datasheet

Rat MMP-24 ELISA Kit (orb567538)

Description	Rat MMP-24 ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	MMP-24
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q99PW6
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat MMP-24(Matrix Metalloproteinase 24) ...