


Product Datasheet

Rat NE ELISA Kit (orb567551)

Description	Rat NE ELISA Kit
Reactivity	Rat
Range	1.563-100ng/ml
Target	NE
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P00774
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Rat
NE (Neutrophil
Elastase)
ELISA Kit