


Product Datasheet

Rat NTXI ELISA Kit (orb567567)

Description	Rat NTXI ELISA Kit
Reactivity	Rat
Range	3.125-200ng/ml
Target	NTXI
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.



Rat
NTXI(Cross
Linked N-
Telopeptide
of T...