


Product Datasheet

Rat OXA ELISA Kit (orb567578)

| | |
|--------------------------|--|
| Description | Rat OXA ELISA Kit |
| Reactivity | Rat |
| Range | 62.5-4000pg/ml |
| Target | OXA |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Application notes | Competitive ELISA, Coated with Antigen |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Sensitivity | 37.5pg/ml |
| Expiration Date | Please enquire. |



Rat
OXA(Orexin
A) ELISA Kit