


Product Datasheet

Rat PDGF-BB ELISA Kit (orb567599)

Description	Rat PDGF-BB ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	PDGF-BB
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q05028
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat PDGF-BB(Platelet-Derived Growth Fact...