


Product Datasheet

Rat PDGF-D ELISA Kit (orb567600)

Description	Rat PDGF-D ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	PDGF-D
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9EQT1
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat PDGF-D (Platelet Derived Growth Facto...