


## Product Datasheet

### Rat PDGFRL ELISA Kit (orb567601)

<b>Description</b>	Rat PDGFRL ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	78.125-5000pg/ml
<b>Target</b>	PDGFRL
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>Q5RJP7</b>
<b>Sensitivity</b>	46.875pg/ml
<b>Expiration Date</b>	Please enquire.



Rat PDGFRL(Platelet Derived Growth Facto...