


Product Datasheet

Rat Pg ELISA Kit (orb567607)

Description	Rat Pg ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	Pg
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	CAS:68-96-2
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
Pg(Progesterone)
ELISA Kit