


Product Datasheet

Rat PTGES ELISA Kit (orb567640)

Description	Rat PTGES ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	PTGES
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	D4AE56
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
PTGES(Prostaglandin
E Synthase, Micr...