


Product Datasheet

Rat ptau ELISA Kit (orb567647)

Description	Rat ptau(phospho Tau Protein) ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	$\rho\tau$
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, cerebrospinal fluid [CSF], Other liquid samples
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
 $\rho\tau$ (phospho
Tau Protein)
ELISA Kit