


Product Datasheet

Rat sgp130 ELISA Kit (orb567667)

Description	Rat sgp130 ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	sgp130
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P40190
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
sgp130(soluble
Glucoprotein
130) ELI...