


## Product Datasheet

### Rat TAP ELISA Kit (orb567689)

<b>Description</b>	Rat TAP ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	156.25-10000pg/ml
<b>Target</b>	TAP
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P00762</b>
<b>Sensitivity</b>	93.75pg/ml
<b>Expiration Date</b>	Please enquire.



Rat  
TAP(Trypsinogen  
Activation  
Peptide) ...