


Product Datasheet

Rat TFF2 ELISA Kit (orb567696)

Description	Rat TFF2 ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/ml
Target	TFF2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q09030
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Rat
TFF2(Trefoil
Factor 2)
ELISA Kit