


Product Datasheet

Rat UCHL1 ELISA Kit (orb567726)

Description	Rat UCHL1 ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/ml
Target	UCHL1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q00981
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Rat
UCHL1(Ubiquitin
Carboxyl
Terminal Hy...