


Product Datasheet

Rat ucOC ELISA Kit (orb567727)

Description	Rat ucOC ELISA Kit
Reactivity	Rat
Range	0.078-5ng/ml
Target	UCOC
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.047ng/ml
Expiration Date	Please enquire.



Rat
ucOC(Undercarboxylated
Osteocalcin) ...