


Product Datasheet

Rat VIP ELISA Kit (orb567733)

Description	Rat VIP ELISA Kit
Reactivity	Rat
Range	7.813-500pg/ml
Target	VIP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01283
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Rat
VIP(Vasoactive
Intestinal
Peptide) E...