


## Product Datasheet

### Rat Caspase-2 ELISA Kit (orb567759)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Rat Caspase-2 ELISA Kit  |
| <b>Reactivity</b>      | Rat  |
| <b>Range</b>           | 0.156-10ng/ml  |
| <b>Target</b>          | Caspase-2  |
| <b>Storage</b>         | 2-8 °C for 6 months  |
| <b>Note</b>            | For research use only  |
| <b>Sample Types</b>    | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| <b>Uniprot ID</b>      | <b>P55215</b>  |
| <b>Sensitivity</b>     | 0.094ng/ml   |
| <b>Expiration Date</b> | Please enquire.  |



Rat  
Caspase-  
2(Caspase-  
2) ELISA Kit