


Product Datasheet

Rat 8-OHdG ELISA Kit (orb567772)

Description	Rat 8-OHdG ELISA Kit
Reactivity	Rat
Range	1.563-100ng/ml
Target	8-OHdG
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, urine, Other liquid samples
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Rat 8-OHdG(8-Hydroxydeoxyguanosine) ELIS...