


Product Datasheet

Rat N-MID-OT ELISA Kit (orb567774)

Description	Rat N-MID-OT ELISA Kit
Reactivity	Rat
Range	1.563-100ng/ml
Target	N-MID-OT
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P04640
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Rat N-MID-OT(N-MID Osteocalcin) ELISA Ki...