


Product Datasheet

Rat E2 ELISA Kit (orb567779)

Description	Rat E2 ELISA Kit
Reactivity	Rat
Range	12.5-800pg/ml
Target	E2
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	7.5pg/ml
Expiration Date	Please enquire.



Rat
E2(Estradiol)
ELISA Kit