


Product Datasheet

Rat 4-HNE ELISA Kit (orb567805)

Description	Rat 4-HNE ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	4-HNE
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	18.75 pg/ml
Expiration Date	Please enquire.



Rat 4-HNE(4-Hydroxynonenal)
ELISA Kit