


Product Datasheet

Rat CXCR4 ELISA Kit (orb567809)

Description	Rat CXCR4 ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	CXCR4
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	008565
Sensitivity	0.039ng/ml
Expiration Date	Please enquire.



Rat
CXCR4(C-X-C
motif
receptor 4)
ELISA ...