


Product Datasheet

Rat GHRL ELISA Kit (orb567821)

Description	Rat GHRL ELISA Kit
Reactivity	Rat
Range	7.8-500pg/ml
Target	GHRL
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9QYH7
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.



Rat
GHRL(Ghrelin)
ELISA Kit