


## Product Datasheet

### Rat Ang1-7 ELISA Kit (orb567828)

<b>Description</b>	Rat Ang1-7 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	7.8-500pg/ml
<b>Target</b>	Ang1-7
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich ELISA, Double Antibody
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	4.688pg/ml
<b>Expiration Date</b>	Please enquire.



Rat Ang1-7 (Angiotensin 1-7) ELISA Kit