


## Product Datasheet

### Rat fT3 ELISA Kit (orb567832)

<b>Description</b>	Rat fT3 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.781-50pg/ml
<b>Target</b>	fT3
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Application notes</b>	Competitive ELISA, Coated with Antigen
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Sensitivity</b>	0.469pg/ml
<b>Expiration Date</b>	Please enquire.



Rat fT3(Free  
Triiodothyronine)  
ELISA Kit