


Product Datasheet

Rat TRPM7 ELISA Kit (orb567841)

Description	Rat TRPM7 ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	TRPM7
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q925B3
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
TRPM7(Transient
Receptor
Potential C...