


## Product Datasheet

### Rat Calprotectin ELISA Kit (orb567843)

<b>Description</b>	Rat Calprotectin ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	6.25-400ng/ml
<b>Target</b>	Calprotectin
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, fecal extracts, Other liquid samples
<b>Sensitivity</b>	3.75ng/ml
<b>Expiration Date</b>	Please enquire.



Rat  
Calprotectin(Calprotectin)  
ELISA Kit