


Product Datasheet

Rat CPP ELISA Kit (orb567863)

Description	Rat CPP ELISA Kit
Reactivity	Rat
Range	3.12-200pg/ml
Target	CPP/Copeptin
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01186
Sensitivity	1.875pg/ml
Expiration Date	Please enquire.



Rat
Copeptin
ELISA Kit