


Product Datasheet

Rat GA ELISA Kit (orb567871)

Description	Rat GA ELISA Kit
Reactivity	Rat
Range	15.6-1000pmol/ml
Target	GA
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	9.375pmol/ml
Expiration Date	Please enquire.



Rat
GA(Glycated
Albumin)
ELISA Kit