


Product Datasheet

Rat LRP1 ELISA Kit (orb567883)

Description	Rat LRP1(LDL receptor-related protein 1) ELISA kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	LRP1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	G3V928
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat LRP1(LDL receptor-related protein 1)...