


Product Datasheet

Rat SFRP4 ELISA Kit (orb567907)

Description	Rat SFRP4 ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	Sfrp4
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9JLS4
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
Sfrp4(Secreted
frizzled-
related prot...