


Product Datasheet

Rabbit ApoA1 ELISA Kit (orb567918)

Description	Rabbit ApoA1 ELISA Kit
Reactivity	Rabbit
Range	31.25-2000ng/ml
Target	ApoA1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P09809
Sensitivity	18.75ng/ml
Expiration Date	Please enquire.



Rabbit
ApoA1(Apolipoprotein
A1) ELISA Ki...