


Product Datasheet

Rabbit ApoE ELISA Kit (orb567923)

Description	Rabbit ApoE ELISA Kit
Reactivity	Rabbit
Range	3.125-200ng/ml
Target	ApoE
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P18287
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.



Rabbit
ApoE(Apolipoprotein
E) ELISA Kit