


Product Datasheet

Rabbit CD40L ELISA Kit (orb567935)

Description	Rabbit CD40L ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/ml
Target	CD40L/TNFSF5
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	G1SKP7
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rabbit
CD40L/TNFSF5(Cluster
of Different...