


## Product Datasheet

### Rabbit CTGF ELISA Kit (orb567943)

<b>Description</b>	Rabbit CTGF ELISA Kit
<b>Reactivity</b>	Rabbit
<b>Range</b>	46.875-3000pg/ml
<b>Target</b>	CTGF
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>G1SIE6</b>
<b>Sensitivity</b>	28.125pg/ml
<b>Expiration Date</b>	Please enquire.



Rabbit  
CTGF(Connective  
Tissue Growth  
Fac...