


Product Datasheet

Rabbit ET-1 ELISA Kit (orb567950)

Description	Rabbit ET-1 ELISA Kit
Reactivity	Rabbit
Range	1.25-80pg/ml
Target	ET-1
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P29560
Sensitivity	0.625pg/ml
Expiration Date	Please enquire.



Rabbit ET-1 (Endothelin 1) ELISA Kit