


Product Datasheet

Rabbit LH ELISA Kit (orb568050)

Description	Rabbit LH ELISA Kit
Reactivity	Rabbit
Range	0.469-30ng/ml
Target	LH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Sandwich ELISA, Double Antibody
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q6IY74
Sensitivity	0.281ng/ml
Expiration Date	Please enquire.



Rabbit
LH(Luteinizing
Hormone)
ELISA Kit