


## Product Datasheet

### Sheep IL-8 ELISA Kit (orb568080)

<b>Description</b>	Sheep IL-8 ELISA Kit
<b>Reactivity</b>	Sheep
<b>Range</b>	15.625-1000pg/ml
<b>Target</b>	IL-8
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P36925</b>
<b>Sensitivity</b>	9.375pg/ml
<b>Expiration Date</b>	Please enquire.



Sheep IL-8(Interleukin 8) ELISA Kit