


## Product Datasheet

### Sheep FSH ELISA Kit (orb568086)

<b>Description</b>	Sheep FSH ELISA Kit
<b>Reactivity</b>	Sheep
<b>Range</b>	3.12-200mIU/ml
<b>Target</b>	FSH
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P01227</b>
<b>Sensitivity</b>	1.875mIU/ml
<b>Expiration Date</b>	Please enquire.



Sheep  
FSH(Follicle-  
Stimulating  
Hormone) ...