


Product Datasheet

Sheep LH ELISA Kit (orb568090)

Description	Sheep LH ELISA Kit
Reactivity	Sheep
Range	0.938-60mIU/ml
Target	LH
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01231
Sensitivity	0.563mIU/ml
Expiration Date	Please enquire.



Sheep
LH(Luteinizing
Hormone)
ELISA Kit