


Product Datasheet

ft3 ELISA Kit (orb568125)

Description	ft3 ELISA Kit
Reactivity	All
Range	0.781-50pg/ml
Target	ft3
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.469pg/ml
Expiration Date	Please enquire.



ft3(Free
Triiodothyronine)
ELISA Kit