


Product Datasheet

FT4 ELISA Kit (orb568126)

Description	ft4 ELISA Kit
Reactivity	All
Range	1.563-100pg/ml
Target	ft4
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.938pg/ml
Expiration Date	Please enquire.



ft4(Free
Thyroxine)
ELISA Kit