


Product Datasheet

Vitamin B9 ELISA Kit (orb568134)

Description	Vitamin B9 ELISA Kit
Reactivity	All
Range	78.125-5000pg/ml
Target	FA/VB9
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



FA/VB9(Folic Acid/Vitamin B9) ELISA Kit