


Product Datasheet

3-NT ELISA Kit (orb568154)

Description	3-NT ELISA Kit
Reactivity	All
Range	1.563-100ng/ml
Target	3-NT
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



3-NT(3-Nitrotyrosine)
ELISA Kit