


Product Datasheet

Cynomolgus Complement C5a Protein (orb570258)

Description	Cynomolgus Complement C5a, Tag Free is expressed from E. coli cells. It contains AA Met 678 - Arg 75...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Complement C5a
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 8.4 kDa. The protein migrates as 11 kDa under reducing (R) condition (SDS-PAGE).
Purity	95%
Source	Cynomolgus Complement C5a, Tag Free (orb570258) is expressed from E. coli cells. It contains AA Met 678 - Arg 751 (Accession # XP_015292262.1).
MW	8.4 kDa
NCBI	XP_015292262.1
Expiration Date	6 months from date of receipt.



Immobilized Cynomolgus Complement C5a, T...